

Small Walking Frame


SET-UP AND OPERATION

⚠ Please read these instructions before use


⚠ Do not use a walking frame to climb stairs or to travel on an escalator

Unfolding / folding the walking frame

Unfold each arm outwards until the lock clicks into place. Check both sides of the frame are firmly locked into position before use.



To fold, press firmly on the single button located in the centre of the cross bar to release the locks.




Then fold each side of the frame inwards until they rest flat against the frame.

⚠ Ideal height

To determine ideal height for your walking frame, stand as straight as possible (preferably wearing an everyday pair of shoes) and look straight ahead with your arms relaxed at your sides.


In this position, the top of the handgrips should sit level with your wrist bone when the walking frame is placed in the correct position around your body. You can measure the distance from your wrist bone to the floor as a helpful reference for height adjustment.



Adjusting the height

Press in the brass locking button to disengage, and turn the leg slightly so the button is no longer aligned with the holes.


Move the leg up or down to your ideal height setting, then turn the leg back to align the locking button with the nearest hole. If unable to achieve an exact height match, make the walking frame slightly higher rather than shorter to avoid stooping.




The locking button must be fully engaged and protruding through the height adjustment hole to ensure a positive lock. Always test that the lock is fully engaged before use.

Using the walking frame

With both hands placed firmly around the handgrips at the top of the walking frame, lift and place the frame to a comfortable arm's length ahead of you, making sure the frame is stable on the ground.



Step in towards the frame, stepping off from your weaker leg (if applicable). Do not step all the way to the front of the walking frame; smaller steps are recommended. Always take small steps when turning.



Regular maintenance

Regularly check the walking frame for damage and wear, paying special attention to the condition of the rubber tips. Replace the rubber tips immediately if you see evidence of rips, tears, cracks or other general wear. If the walking frame is damaged in any way, stop using the walking frame immediately and contact your authorised dealer or the manufacturer.

Care & cleaning

To clean the walking frame, wipe with a mild soap and warm water solution and rinse well. Dry thoroughly with a clean cloth.

Avoid harsh, abrasive cleaning agents or utensils (such as a wire brush) as these may damage the anodised finish, or plastic and rubber components.

Warranty

Evadale Healthcare Pty Ltd warrants this product for a period of 12 months. Please refer to our Warranty Statement for full details. Please contact your authorised dealer or Evadale Healthcare for a list of available spare parts.

Specifications

Product Code	ED1461
Maximum User Weight	110 kg
Product Weight	2.5 kg
Overall Width	600 > 670 mm
Overall Height	787 > 965 mm
Overall Depth	460 > 500 mm
Depth When Folded	100 mm
Width Between Handles	430 mm
Warranty	1 year
ARTG	176830

Made in China



www.evadale.com.au

**For more information
or to order, freecall
1300 883 194**

©2012 Evadale Healthcare. All rights reserved.