


Small Walking Frame



Compact folding



A durable yet lightweight aluminium walking frame that is ideal for everyday use. It features height adjustable legs with replaceable rubber tips and foam grip handles for added user comfort.



The frame folds flat to 100 mm for easy transport and storage via a single button folding mechanism.

OVERALL WIDTH	600 > 670 mm
OVERALL HEIGHT	787 > 965 mm
OVERALL DEPTH	460 > 500 mm
DEPTH WHEN FOLDED	100 mm
BETWEEN HANDLES	430 mm
WEIGHT	2.5 kg

CODE	ED1461
RAP	AP12
MATERIAL	Aluminium
ARTG	176830