


The ROHO[®] Insert and Leveling Pads

Insert Pad


Foam insert pad(s) are available in both standard and custom configurations and may be snapped to ROHO[®] DRY FLOATATION[®] Mattress Section(s) to provide patients a localized approach to wound management when cost restraints are required.

MODEL	DESCRIPTION
IP2134	(21" x 34") Insert Pad

Customized foam configurations are available.



Insert pad shown with three ROHO mattress sections, positioned to support the head.




Insert pad shown with two ROHO mattress sections, positioned to support the head and foot area.

NOTE: Foam inserts are designed to be used only in areas of low risk. For high risk clients, we recommend the use of the complete ROHO DRY FLOATATION Mattress System.

Leveling Pad

For protection in the ischial area, the leveling pad allows the mattress section to fit snugly into a cut-out center for no-slip protection.

MODEL	DESCRIPTION
LP3480	(34" x 80") Leveling Pad



One mattress section, which is centered in a precut plastic Staphcheck[®] covered foam leveling pad. This configuration protects the mid-torso area.

ROHO leveling pads and inserts are constructed of materials that voluntarily meet Federal and State of California specifications for flammability. Each pad is covered with a special Staphcheck[®] material that is both *antibacterial* and *flame retardant*. The material is easily wiped clean and is very durable.



the ROHO group
shape fitting technology[®]

100 North Florida Avenue • Belleville, IL 62221-5429 • USA • www.therohogroup.com
U.S. & Canada 1-800-851-3449 • 618-277-9150 • Fax 618-277-6518
Outside the U.S. & Canada • 1-618-277-9173 • Fax 1-618-277-6518