

# Contoured Wheelchair Cushion

Gently contoured to fit  
the bottom and thighs

Enhanced comfort and  
maximum pressure redistribution



Combining a high density foam support base, with a super soft memory foam top layer, this cushion is gently contoured to fit the bottom and thighs, offering maximum pressure redistribution and enhanced comfort to minimise the risk of pressure injury.

This combination provides the benefits of a stable sitting surface and even pressure redistribution, without constricting blood circulation, lessening the risk of pressure injury.

This cushion comes with a durable waterproof cover that has a zipper along the rear edge for easy removal. The base of the cover is made from non-slip

material that helps to stop the cushion from sliding forward on the chair. The cover includes a skirted flap that lays over the zipper to increase infection control – an extra flap of material that helps to cover the zipper to stop fluid seeping into the layers of the cushion.

The cover is machine washable for easy cleaning and hygiene.

<b>CODE</b>	BL7265
<b>HEIGHT</b>	105 > 65 mm
<b>WIDTH</b>	430 mm
<b>DEPTH</b>	410 mm
<b>CAPACITY</b>	130 kg
<b>WARRANTY</b>	1 year

