

Coraille

Shell chair for indoor / outdoor usage

total  mobility

2 big rear wheels and 2 smaller front wheels.

Aluminium frame.

One-piece spring mounted footplate of high resistance ABS synthetic material.

Pneumatically adjustable chair angle and leg supports

5-point safety harness

Padded headrest

Footbrake at the back

Seat height adjustable from 45cm - 56cm




Footbrake on rear wheel














OPTIONAL removable/foldable table
(telescopic adjustment 38 - 48cm)



OPTIONAL B52 Serum holder



OPTIONAL Separate footrests

|  |  |  |  |  |  |  |  |  |  |  |
|---|---|---|---|---|---|---|---|---|---|---|
| Width | Seat Width | Armrest Height | Seat Height | Total Height | Backrest Height | Seat Incline | Seat Depth | Total Length | Unit Weight | Weight Capacity |
| 67cm | 38cm - 48cm | 23cm - 27cm | 45cm - 56cm | 119cm - 124cm | 75cm - 79cm | -4° - 36° | 46cm | 112cm | 36kg | 120kg |
| 790 (XXL) | 55cm - 60cm | 23cm | 52cm - 56cm | 124cm | 78cm | | 50cm | 112cm | 37kg | 150kg |