


THE MERIT MATTRESS

Our original design, and still firm favourite.

The Merit is ideal for low care applications where support and comfort is a must. Often used in bespoke situations including high humidity conditions and custom size homecare beds, the Merit will always ensure quality rest and care. The Merit is also an excellent choice for short term hire applications where durability is a must.

STANDARD CONFIGURATIONS INCLUDE


- In-built convolute airflow design for pressure care and comfort
- Single core, deluxe grade foam — often used to support Air Mattress systems



FEATURES

DELUXE CONVOLUTE FOAM

- Provides excellent durability
- Unique deep profile assists with internal air circulation and temperature management
- Supportive and comfortable



WATER-BASED EVAPORATIVE LAMINATION

- Completely odour free
- Environmentally friendly
- Incredible Adhesion



OPTIONS

1 PROFILE STABILISING SIDES (THIS OPTION SHOWN ON MAIN IMAGE)

- Enables ease of primary patient care
- Facilitates patient transfers
- Unique hinge points allow mattress to easily conform to varied bed positions
- Encourages central patient positioning



2 SAFETY SIDES

- Provide unobtrusive support
- Help prevent high falls risk patients rolling out of bed
- Gently encourages central patient positioning to maximise pressure care

SPECIFICATIONS

The Merit can be adapted to suit any bed size, configuration or profiling function.

STANDARD SIZING	length: 1980mm
	2030mm
thickness:	50mm
	100mm
	(single)
widths:	880mm
	1050mm (king single)
	1150mm
	1350mm

WEIGHT RANGE*	Single Core 100-150kg
	Covolute 130-220kg (depending on selected width)

WARRANTY 5 years**

*Optimal patient weight range for pressure care and comfort.

**Warranty refers to any manufacturing or product defects. Due to its nature, foam will soften and its resilience will decrease with usage and age. Contact Forté Healthcare Solutions for further specific information.

PRESSURE CARE COVER

FEATURES

1 MULTI-STRETCH PRESSURE CARE COVER UPPER

- Minimises shear
- Anti-Microbial and resistant to body fluids
- Facilitates Immersion of pressure points

2 HEAVY DUTY PVC BASE

- Heavy duty protection for base of cover
- Minimises shear forces during transit by sliding easily on bed base
- Anti-microbial and resistant to body fluids

3 INTEGRATED ZIP WITH PROTECTION FLAP

- Allows removal of cover for cleaning
- Flap prevents ingress of fluid through zip
- Fully welded construction

COVER OPTIONS

- Full multi-stretch pressure care fabric
- High frequency sealed and welded seams
- Fully enclosed and welded cover design



MAKING PATIENT CARE EASY



P. 1300 73 11 63
W. fortehealthcare.com.au
E. sales@fortehealthcare.com.au