

GEC5 - ECONOMY



The GEC5 enables one person to evacuate a physically impaired person safely and easily down stairs in the event of an emergency or when lifts can not be used.

Comes complete with operating instructions that are visible to the operator at all times, a free dust cover and wall mounting brackets.

The GEC5 is our economy model with many of the features of the rest of the GLOBEX range.

Operation:	Down stairs and wheelchair transportation
Construction:	Aluminium
Open height (mm):	1600
Folded height (mm):	1040
Width (mm):	500
Folded dimensions (mm):	210
Weight (kg):	10.0
Loading weight (kg):	150



Fold away secondary

Friction controlled high resistance belts



Robust easy grip handle

Adjustable headrest /back support with operator instructions

Telescopic handle with automatic securing non pinch clips