



oscar
FURNITURE

M5 Bariatric Lift-Recline Chair

Style, performance and comfort



Features:

- 320kg seat lift capacity
- 75kg footrest lift capacity
- Pressure care memory foam as standard
- Triple-action motor control to allow for more comfortable chair positions
- "Made to measure" options available

Options & Modifications:



Specifications:

Size	Seat Depth	Back Height	Arm Height	Footrest	Seat Width	Seat Height
650	530mm	690mm	190mm	950mm	650mm	500mm
650/70	600mm	770mm	190mm	1020mm	650mm	570mm
850	530mm	690mm	190mm	950mm	850mm	500mm
850/70	600mm	770mm	190mm	1020mm	850mm	570mm