



## Nordic Mobicline

### Description:

A versatile mobile chair that also reclines. This chair features a fully welded tubular frame and incorporates a large activity tray and slide out footrest.

### Options available:

150mm castors  
Directional locking castors  
Wingback head supports

### Fabric Required

4.5 metres

### Dimensions:

Overall Width: 700mm  
Seat Depth: 490mm  
Seat Height (from footrest): 445mm  
Seat to top of back: 650mm  
Seat to top of arm: 190mm  
Length (extended): 1410mm

### Lead Time:

Allow 4 weeks in standard fabric

### Warranty:

Metal frame: 5years

### Points of Note:

Mobile recliner chair  
Fully welded rigid frame  
Lockable 100mm castors  
Large removable activity tray with holder  
Slide out footrest plate  
3 segment back cushion  
Removable cushions for ease of cleaning

Updated April 2009