

## Ottobock - Start M4 XXL



### Stability

We designed the very stable **Start M4 XXL** especially for large, strong persons. It will support a load of up to 160 kilograms.

The frame has a sturdy double cross brace that can withstand extreme loads. In addition, the back tubes with lumbar angle are reinforced.

The upholstery is also extremely stress-resistant.

Attributes	Details
<b>Supplier</b>	<b>Total Mobility</b>
<b>Classification</b>	Medium Active
<b>Propulsion</b>	Self-Propelling
<b>Max User weight</b>	<b>160 Kg</b>
<b>Seat Width (cm)</b>	45.5 - 58
<b>Seat Depth (cm)</b>	43 - 49
<b>Seat Height (cm)</b>	<b>39.5 - 52 (Front); 37 - 52 (Rear)</b>
<b>Seat Type</b>	Sling / Seat Inclination of 3°
<b>Backrest Type</b>	Sling / Angle Adj
<b>Backrest Height (cm)</b>	40 - 47.5
<b>Backrest Material</b>	Durable and Resistant Upholstery
<b>Push Handles</b>	Fixed
<b>Legrest</b>	Calf Supports / Elevating
<b>Footrest</b>	Flip Up / Swing Away / Angle Adj
<b>Armrest</b>	Detachable / Desk Side Panels / Height Adj
<b>Available Colours</b>	3 Frame Colours Available
<b>Tyres</b>	Solid / Pneumatic
<b>Safety Aspects</b>	Lap-Strap / Anti-tip Bars / Optimized Wheel Lock
<b>Wheelchair Weight</b>	From 17.5 Kg
<b>Other</b>	Integrated Tip Assist / Solid Double-cross Brace

Unit 10/1-11 Burns Road,  
Heathcote, NSW, 2233, Australia.  
Sales Office-West Pennant Hills



Tel: (02) 9520 1866

Fax: (02) 9520 4345

Email:

[info@totalmobility.com.au](mailto:info@totalmobility.com.au)

[www.totalmobility.com.au](http://www.totalmobility.com.au)

ABN: 45 122 616