



SOLID WHEELCHAIR RAMPS

“Full Platform provides excellent stability

“Durable welded aluminium fabrication

“Anti slip traction surface

“Non slip surface

“Full 76cm or 91.5 cm wide platform provides excellent stability for both 3 & 4 wheel model scooters and wheelchairs

“272kg carrying capacity

“Available in 91cm (3ft) 122cm (4ft) and 152cm (5ft) lengths and 2 widths of 76cm or 91.5cm

“Weight carrying Capacity 136kg on one axel or 272kg on 2 axels