



EVACUATION CHAIRS  
AUSTRALIA

TEL: +61 3 9001 5806

1300 669 730

FAX: +61 3 9011 9693

EMAIL: [office@evacuationchairs.com.au](mailto:office@evacuationchairs.com.au)

WEB: [www.evacuationchairs.com.au](http://www.evacuationchairs.com.au)

## Spencer Pro-Skid Evacuation Chair



The Pro-Skid is meant for use in public buildings because, when necessary, it is ready to be used without any compromise, and not only in the case of evacuation.

Designed to maximize the manoeuvrability in cramped spaces and on stairs and to provide easy transport of personnel to and from stairs and towards the outside of the building.

On the Pro-Skid, the back rest and the foot rest extend to become support handles with an ergonomic technique of lifting, thus improving the view of the stairs and feet for the operator, when obstacles on the stairs necessitate the lifting of the chair.

### More information



### Standard features include:

- The chair has a caterpillar belt system that allows movement down the stairs with the help of one operator only.
- The frame's tubes have a minimum diameter of 25 mm up to a maximum of 35 mm in order to improve stability, strength and durability.
- The front wheels have a brake system which blocks both horizontal and vertical rotation to stop movement, even on steep slopes.
- The chair has a headrest with head-fixing belt system in order to improve safety.
- Both headrest and seat are detachable, padded and made out of vinyl material coating 600D, antibacterial and completely filled in order to improve comfort, hygiene and durability.
- The chair has posterior integrated carrying handles that attach at the support frame and to the hub between the chassis and the support frame, anterior telescopic handles to make possible also the use of the chair to ascend a stair (2 person operation), say from a basement.
- Integrated instruction provided on the head rest.

### Technical Data

Weight	18Kg
Maximum load	150Kg
<b>Chair folded</b>	
Height	130cm
Width	55cm
Depth	45cm
<b>Chair fopen</b>	
Height	160cm
Width	55cm
Depth	80cm