

# Custom Sling Design Script

Date \_\_\_\_\_ Therapist \_\_\_\_\_ Client \_\_\_\_\_

## 1- CHOOSE A BASE SLING DESIGN (recommended) plus your preferred fixing.

Base Sling Design (loop fixing only)		Base Sling Design	Loop	Clip	E-Clip	Dual (Clip)	Dual (E-Clip)			
General Purpose		Pivot								
General Purpose Head Support		Pivot Care								
Care		Standing Transfer								
Care Head Support		Super Soft								
Junior		Other								
Amputee		<b>Base Material Type</b>								
Deluxe Leg							Closed Fabric		Open Mesh	
Walking										

## 2- CHOOSE A BASE SLING SIZE. If you have selected "Other" please provide generic dimensions (A-D) plus any other relevant dimensions in the sketch box provided on the next page

Base Sling Size						
Extra Small		Small		Medium		Large
X Large		XX Large		Other		

## 3- CHOOSE A STANDARD MODIFICATION.


### MOD 001 - add/change torso belt

For standard Torso Belt lengths/widths, please see page 6 of sling brochure.

Torso Belt					
Buckle		Neoprene/Buckle Over		Custom Length	
Velcro		Neoprene/Velcro		Custom Width	
Comment					


### MOD 002 - add/reduce head support

Head Support	
Add Head Support	cm
Reduce Head Support	cm




### MOD 003 - add/reduce body length

Body Length	
Add Body Length	cm
Reduce Body Length	cm



### MOD 004 - add/reduce body width

Body Width	
Add Body Width	cm
Reduce Body Width	cm



NB - Body width adjustment will not affect generic external profile of sling.

