


Bariatric Electric Lift & Recliner Chairs

ndis Official Provider



140 - 250 KG Capacity



**SIT, RECLINE,
REST AND
STAND**

**ELEVATED FOOTREST
85mm - 200mm
ABOVE SEAT HEIGHT**

**STEEL SEAT
FRAME WITH
STEEL SPRINGS**

 **ambassador**





Ambassador is a market leader in the manufacture and supply of superior lift & recline chairs for the healthcare market.



We've utilised our expertise in this field to blend exceptional comfort and advanced technologies, designing the ultimate bariatric recliner chair with exclusive mechanisms and components that are extra strong and durable to provide greater lifting capacity than any other lift and recliner chair in this category.

Ambassador chairs enable more independence in the home and numerous benefits for the caregiver or nursing staff.

UNIQUE MECHANICS AND STEEL FRAME

- Provides a strong foundation for each chair
- Using steel instead of timber or chipboard, this guarantees the chair lasts longer without sagging or losing its shape
- The use of bolts improves durability and enables easier repairs
- Longer life span allows for a longer industry warranty
- Flexibility for versatility and improved individual comfort

UNIQUE MOTOR & ELECTRICS FOR ADJUSTABLE POSITIONS

- Fine-tune seating for optimum comfort and easy adjustment
- Results in improved posture and comfort
- Recline, sit up or stand at the touch of a simple button
- Provides simultaneous movement of backrest and footrest with no manual adjustments necessary
- Reliable and strong performance for reassurance
- Alleviate the reliance and risk of injury for the caregiver or nursing staff.
- Medical grade LINAK electronics suitable for patients with pace makers.

FEATURES

Removable chair covers, seat and backrest cushions for ease of cleaning or replacement

Headrest cushion for comfort

Choice of 2-3
lumber supports

Single hand control to recline, sit or stand up at the touch of a simple button

Elevated footrest above seat height for circulation and recovery



Bronze bushes in all lifting points

Upholstered in the highest quality soft to touch fabric. In a choice of 4 standard colours and optional extras.

Steel seat frame with NO SAG springs for maximum strength, support and comfort



Bolted steel frame improves durability and longevity and enables easier repairs

Rechargeable back up battery operating system

European Medical Grade Electronics
– Pace maker compliant



DIMENSIONS

	Leg & back configuration	Leg Rest mm	Lift Capacity KG	Seat height cm	Seat depth cm	Seat width cm	Arm height cm	Back height cm	Overall width cm	Overall height cm	Overall length cm	Back Cushions
Premier A1	Simultaneous leg and back	200 mm higher than the body	140	50	48	50	17	77	80	115	95	2
Premier A2	Independent back and leg	85mm higher than the body	160	50	48	50	17	77	80	115	95	2
Premier A3	Independent back and leg	85 mm higher than the body	250	50	55	68	20	78	102	115	100	2
Senator A1	Simultaneous leg and back	200 mm higher the body	140	50	48	50	17	77	80	115	95	3

STANDARD COLOURS



All models come in a selection of 4 standard colours.

OPTIONAL EXTRAS

Additional fabrics, custom sizes, massage and lumber options available at additional cost and lead-time.

SENATOR A1

PREMIER A1/A2

PREMIER A3

WIDE VARIETY
OF PREMIUM
SHANN LEATHER
OPTIONS
AVAILABLE

FOR YOUR NEAREST RETAILER AND MORE INFORMATION:

Free Call: 1800 807 896 Tel: 02 9622 1882

www.ambassadorproducts.com.au