

küschall®  
designed for life















*Simply extreme and superlight.*

The *K*<sub>SL</sub>

# The K<sub>SL</sub>




## Technical Data

 Seat width: 340-420 mm in increments of 20 mm	 Backrest angle: 74°/78°/82°/86°/90°
 Seat depth: 375 mm - 475 mm in increments of 25 mm	 Total width: SW + 170 mm
 Seat height: front: 450 - 510 mm rear: 400 - 480 mm in increments of 10 mm	 Total length: 75°: approx. 800 mm 90°: approx. 730 mm
 Backrest height: 270 - 390 mm in increments of 15 mm	 <b>Total weight unloaded: approx. 6,7 kg</b>
 Knee-to-heel length: 400-480 mm in increments of 10 mm	 Max. user weight: 100 kg
 Center of gravity: 80 mm/ 95 mm/ 110 mm/ 125mm/ 140 mm	 <b>Transport weight w/out rear wheels, Approx. 4,5 kg</b>

All measurements and weights indicated are based on seat width 400 mm and seat depth 400 mm in the lightest and most compact configuration possible. Weight and measurements can change with the different configurations.

## Options

Made in Switzerland 




These are only a few of the many options available for this model. Contact your nearest küschall specialist for more information.

## Colours



## CPS (Colour Power Set)



The Colour Power Set (CPS) applies on the following parts and options: High Performance rear wheel hubs with quick release protection, fork, Starec front wheels and castor housing caps.

# A Masterpiece; Lightness and design in perfect harmony.

The lightest chair on the market: 4.5 kg (without rear wheels)



*Freedom through supreme lightness.*

**Oliver Anthofer,**

*Freedom is key. The KSL by küschall is the ultimate culmination of design, performance and extreme lightness giving me the perfect means for complete independence. With it I can achieve almost anything I set my mind to.*

*Life is a succession  
of moments. To live  
each one is to  
succeed.* Corita Kent